



# 一次性使用旋转型活体取样钳

DISPOSABLE ROTATABLE BIOPSY FORCEPS

全国**首创** 单独**取证**




## 反复开合 同步**旋转**

- 反复开合, 360°同步旋转 ✓
- 任意角度取样, 精准抓取 ✓
- 内轴旋转, 减少镜子损耗 ✓
- 行业一流制作工艺 ✓



A款



B款

锐意进取 敬天爱人

江苏锐天医疗科技有限公司 Jiangsu Ruitian Medical Technology Co., Ltd



## ·产品特点

- 1、反复开合，360°同步旋转，灵活应对病灶部位在任意角度取样，抓取精准，操作便捷。
- 2、内轴旋转，减少镜子损耗。
- 3、活检钳切口锋利度，吻合度，咬合力，操作性，制造工艺都已达到行业一流水平。

## ·专家评价

- 1.必要性：一次性使用旋转型活体取样钳的钳头可360°旋转，针对取活检角度比较困难、位置不是特别好的时候，可通过旋转钳头角度来解决这一问题，抓取更加精准。
- 2.可行性：通过调整角度来抓取活检，操作更加简单方便，可大幅减少手术操作时间，提高活检效率。
- 3.安全保障：本产品经过各大医院临床试用，深受主任医生欢迎。产品工艺成熟，安全可靠，四连杆结构抓取更有力。钳头有独特的切口设计，极大减少了取样时的出血量，各方面性能均优于普通活检钳。

规格型号 Model Number	直径(mm) Diameter	长度(mm) Length	钳道(mm) Forceps channel	特征 Features
QFB-A-23-1600	2.3	1600	≥2.8	胃镜、不包塑、不带针 Gastroscope, uncoated, without spike
QFB-A-23-1900	2.3	1900	≥2.8	胃肠镜、不包塑、不带针 Gastroscope, uncoated, spike
QFB-A-23-2300	2.3	2300	≥2.8	肠镜、不包塑、不带针 Endoscope, uncoated, without spike
QFB-C-23-1600	2.3	1600	≥2.8	胃镜、包塑、不带针 Gastroscope, coated, spike
QFB-C-23-1900	2.3	1900	≥2.8	胃肠镜、包塑、不带针 Endoscope, uncoated, without spike
QFB-C-23-2300	2.3	2300	≥2.8	肠镜、包塑、带针 Endoscope, coated, spike
QFB-A-18-1000	1.8	1000	≥2.0	支气管镜、不包塑、不带针 Bronchoscope, uncoated, without spike
QFB-A-18-1200	1.8	1200	≥2.0	支气管镜、不包塑、不带针 Bronchoscope, uncoated, without spike
QFB-A-18-1600	1.8	1600	≥2.0	超细胃镜、不包塑、不带针 Superfine gastroscop, uncoated, without spike
QFB-C-18-1000	1.8	1000	≥2.0	支气管镜、包塑、不带针 Bronchoscope, coated, without spike
QFB-C-18-1200	1.8	1200	≥2.0	支气管镜、包塑、不带针 Bronchoscope, coated, without spike
QFB-C-18-1600	1.8	1600	≥2.0	超细胃镜、包塑、不带针 Superfine gastroscop, coated, without spike



地 址：江苏常州湖塘科技产业园B3幢  
电 话：(0519) 8366 8168 8356 8168  
136 0650 9326 180 0150 5866  
传 真：(0519) 8376 8168 邮 编：213161  
邮 箱：gm@realsky-medical.com  
网 址：http://www.realsky-medical.com

Address: B3 Building, Science and Technology Industrial Park  
Hutagn Changzhou, Jiangsu  
Tephone: (0519) 8366 8168 8356 8168 136 0650 9326 180 0150 5866  
Fax: (0519) 8376 8168 Zip code: 213161  
E-Mail: gm@realsky-medical.com  
Web: http://www.realsky-medical.com